

Digirolebs

1.3

DICOM Conformance Statement

V 1.0

© medigration GmbH
All rights reserved

medigration GmbH, Schuhstr. 30, 91052 Erlangen
Tel: +49 (0) 9131/69087-40 , Fax: +49 (0) 9131/69087-50,
www.medigration.de

Table of Contents

| | |
|---|----------|
| 1 Introduction | 5 |
| 1.1 Purpose | 5 |
| 1.2 Scope..... | 5 |
| 1.3 Definitions, Acronyms and Abbreviations..... | 6 |
| 1.3.1 Definitions | 6 |
| 1.3.2 Acronyms and Abbreviations | 6 |
| 1.4 References..... | 6 |
| 2 Implementation Model | 7 |
| 2.1 Application Data Flow Diagram..... | 7 |
| 2.2 Functional Definition of Application Entities | 8 |
| 2.3 Sequencing of Real World Activities | 8 |
| 3 Digiroebs Application Entity Specification | 9 |
| 3.1 Association Establishment Policies..... | 9 |
| 3.1.1 General | 9 |
| 3.1.2 Number of Associations | 9 |
| 3.1.3 Implementation Identifying Information | 9 |
| 3.2 Association Initiation Policy..... | 9 |
| 3.2.1 Operator Initiated Communication Verification..... | 10 |
| 3.2.1.1 Associated Real-World Activity: User Echo (Dicom Option)..... | 10 |
| 3.2.1.2 Proposed Presentation Contexts | 10 |
| 3.2.1.3 SOP Specific Conformance | 10 |
| 3.2.1.4 Association Termination..... | 10 |
| 3.2.2 User Initiated Image Send..... | 11 |
| 3.2.2.1 Associated Real-World Activity: User Send..... | 11 |
| 3.2.2.2 Proposed Presentation Contexts | 11 |
| 3.2.2.3 SOP Specific Conformance for DX Storage SOP Class | 11 |
| 3.2.2.3.1 Status codes for C-STORE | 11 |
| 3.2.2.3.2 Supported IOD modules | 12 |
| 3.2.2.3.3 Supported DICOM Elements..... | 13 |
| 3.2.2.4 Association Termination..... | 20 |

| | |
|--|-----------|
| 3.2.3 User initiated Print (Dicom Option)..... | 21 |
| 3.2.3.1 Associated Real-World Activity..... | 21 |
| 3.2.3.2 Proposed Presentation Contexts | 21 |
| 3.2.3.3 SOP Specific Conformance | 21 |
| 3.2.3.3.1 Basic Grayscale Print Management..... | 21 |
| 3.2.3.3.2 Basic Film Session..... | 22 |
| 3.2.3.3.3 Basic Film Box | 22 |
| 3.2.3.3.4 Basic Grayscale Image Box..... | 24 |
| 3.2.3.3.5 Printer | 25 |
| 3.2.4 Operator initiated query for Modality Worklist | 26 |
| 3.2.4.1 Associated Real-World Activity..... | 26 |
| 3.2.4.2 Proposed Presentation Contexts | 26 |
| 3.2.4.3 SOP Specific Conformance | 27 |
| 3.2.4.3.1 Status codes for C-FIND | 27 |
| 3.2.4.3.2 Supported DICOM Elements..... | 27 |
| 3.2.5 Operator initiated Sending of Study Results | 28 |
| 3.2.5.1 Associated Real-World Activity..... | 28 |
| 3.2.5.2 Proposed Presentation Contexts | 28 |
| 3.2.5.3 SOP Specific Conformance | 29 |
| 3.2.5.3.1 N-CREATE attributes..... | 29 |
| 3.2.5.3.2 N-SET attributes | 30 |
| 3.3 Association Acceptance Policy | 33 |
| 4 Communication Profiles | 34 |
| 4.1 Supported Communication Stacks..... | 34 |
| 4.1.1 TCP/IP Stack | 34 |
| 4.1.1.1 Physical Network Media Support..... | 34 |
| 5 Configuration | 35 |
| 6 Support of Extended Character Sets | 36 |

List of Figures

Figure 1: DigiRoebis Implementation Model7

List of Tables

Table 1: Supported DICOM SOP Classes and Roles9

Table 2: Proposed presentation context for verification10

Table 3: Proposed presentation contexts for operator initiated image send11

Table 4: Behavior when receiving C-STORE response status codes (operator initiated)12

Table 5: IOD modules used by DigiRoebis.....13

Table 6: Proposed presentation context for print21

Table 7: SOP Classes of Basic Grayscale Print Management Meta SOP Class.....21

Table 8: Supported DIMSE Services for Basic Film Session22

Table 9: Supported SOP Class Elements for Basic Film Session22

Table 10: Supported DIMSE Services for Basic Film Box23

Table 11: Supported SOP Class Elements for Basic Film Box24

Table 12: Supported DIMSE Services for Basic Grayscale Image Box24

Table 13: Supported SOP Class Elements for Basic Grayscale Image Box25

Table 14: Supported DIMSE Services for Printer26

Table 15: Supported SOP Class Elements for Printer26

Table 16: Proposed presentation contexts for operator initiated query for modality worklist.....27

Table 17: Behavior when receiving C-STORE response status codes (operator initiated)27

Table 18: Dicom attributes for Modality Worklist28

Table 19: Proposed presentation contexts for operator initiated sending of study results29

Table 20: DICOM attributes sent in a N-CREATE-RQ30

Table 21: DICOM attributes sent in a N-CREATE-RQ33

1 Introduction

1.1 Purpose

A DICOM Conformance Statement is intended to describe which components, optional components or extensions of the DICOM standard are supported by a particular implementation. The Conformance Statement of one implementation can be compared with the Conformance Statement from another implementation to determine which capabilities are commonly supported.

DICOM does not, by itself, guarantee interoperability. Furthermore, the identification of common capabilities by comparing DICOM Conformance Statements is also not sufficient to guarantee connectivity between two devices.

A DICOM Conformance Statement cannot replace validation and cross-vendor testing with other devices. Validation and cross-vendor testing are still required to ensure that both devices are performing as intended.

The reader should be aware of a number of important issues:

- Even when comparing this Conformance Statement with the Conformance Statement of another device indicates that connectivity is possible, the system integrator is responsible for carrying out test procedures to ensure that the required connectivity is actually met.
- Neither the DICOM Standard nor this Conformance Statement can ensure interoperability when integrating devices from different vendors. It is the system integrator's responsibility to ensure that the application requirements of all devices within the complete system are met.
- The DICOM standard undergoes continual review and improvement in order to meet changing requirements. Corrections, extensions and additional services are added from time to time. Medigration reserves the right to make changes to the product described in this conformance statement in order to cover changes in the DICOM standard. Readers should be aware that connected devices should also follow changes in the DICOM standard in order to retain connectivity.

The intended audience for this Conformance Statement is hospital technical staff, system integrators and software engineers. The reader is assumed to have good understanding of the DICOM standard.

1.2 Scope

This conformance statement describes the DICOM capabilities of the medigration product **Digirolebs**. The **Digirolebs** is a digital radiography system producing DICOM images. It is specifically designed to be integrated into a PACS environment containing other DICOM storage devices like archives and workstations. It supports those DICOM services needed to store its images for archival and to print them out

on film or paper.

1.3 Definitions, Acronyms and Abbreviations

1.3.1 Definitions

| | |
|-------------------|---|
| System Integrator | A person or organization responsible for integrating devices into a new or existing system. The System Integrator takes responsibility for ensuring that the system works as a whole. |
|-------------------|---|

Other definitions can be found within the different parts of the DICOM standard [1].

1.3.2 Acronyms and Abbreviations

| | |
|------|--|
| AE | A pplication E ntity |
| DCO | D ICOM C omposite O bject. A DICOM object such as an image, overlay, lookup-table, waveform, presentation state or radiotherapy plan which can be stored using the Storage Service Class. |
| DR | D igital R adiography |
| DX | D igital X -Ray |
| MPPS | M odality P erformed P rocedure S tep |
| MWL | M odality W orklist |
| SCU | S ervice C lass U ser |
| SW | S oftware |

Other acronyms and abbreviations used within this document are defined within the different parts of the DICOM standard [1].

1.4 References

- [1] DICOM, PS3.(1-18)-2004, National Electrical Manufacturers Association, 1300 N. 17th Street Rosslyn, Virginia 22209, USA.

2 Implementation Model

DigiRoebis is a device, which produces Digital X-Ray (DX) images. These images are stored on local magnetic disks and can be sent over a network interface for storage or can be printed on a hardcopy device. DigiRoebis maintains a database of summary information about the locally stored images.

DigiRoebis also acts as a Modality Performed Procedure Step SCU to send study results each time after an acquisition is performed and can request Modality Worklists.

2.1 Application Data Flow Diagram




Figure 1: DigiRoebis Implementation Model, *Dicom Options are shown in dashed line

Figure 1 illustrates the relationships between the DigiRoebis Application Entity (AE) and its associated Real-World Activities. The **Remote Real-World Activities** are shown on the right and the **Local Real-World Activities** are shown on the left.

Send to Archive is an activity performed by DigiRoebis to send images to a DICOM device to be stored by the **Storage** remote activity.

Query is an activity performed by Digiroebs to receive patient lists from a remote device provided by the **Query Worklist** activity.

User Echo is an activity performed by Digiroebs to verify the ability of a remote device to respond to DICOM messages. Echo messages will be sent upon operation request. They are responded by the **Echo Provider** activity on the remote device.

Send study result is an activity performed by **Digiroebs** to send Modality Performed Procedure Step notifications to a remote device. These notifications are stored by the **Receive study result** activity.

2.2 Functional Definition of Application Entities

The Digiroebs software acts as a single Application Entity (AE), which is able to send images for storage, study result notifications and echo requests and to query a remote worklist database for patient lists.

Digiroebs acts as an SCU of the following DICOM Service Classes:

- Verification
- Storage
- Study Management
- Basic Worklist Management
- Print Management

2.3 Sequencing of Real World Activities

No sequencing of Real-World activities are relevant.

3 Digiroebs Application Entity Specification

Digiroebs provides standard conformance to the following Service Classes by supporting the SOP Classes and roles listed in Table 1.

| SOP Class Name | UID | Role |
|--|-----------------------------|------|
| Verification | 1.2.840.10008.1.1 | SCU |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | SCU |
| Digital X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | SCU |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | SCU |
| Basic Grayscale Print Management | 1.2.840.10008.5.1.1.9 | SCU |

Table 1: Supported DICOM SOP Classes and Roles

3.1 Association Establishment Policies

3.1.1 General

All relevant DICOM communication parameters (AE Titles, hostnames or IP addresses, port numbers, etc.) are configurable. See section 5 for more information on configurable parameters. A maximum PDU size of 16KB will be offered when establishing associations.

3.1.2 Number of Associations

Digiroebs does not accept any DICOM associations.

3.1.3 Implementation Identifying Information

| | |
|------------------------------|-----------------------------|
| Implementation Class UID: | 1.2.276.0.7230010.3.0.3.5.4 |
| Implementation Version Name: | OFFIS_DCMTK_354 |

3.2 Association Initiation Policy

Digiroebs will initiate associations in the following situations:

- When instructed by an operator (via the user interface) to verify communication. (Dicom Option)
- When instructed by an operator (via the user interface) to send images to a remove device.
- When instructed by an operator (via the user interface) to print images to a remote device. (Dicom Option)
- When instructed by an operator to send study results to a remove device.
- When instructed by an operator to query for modality worklists.

3.2.1 Operator Initiated Communication Verification

3.2.1.1 Associated Real-World Activity: User Echo (Dicom Option)

An operator can - via a graphical user interface - initiate a test to verify communication. The associated local real-world activity is **User Echo** and the remote real world activity is **Echo Provider**. The communication verification test is considered successful if an association can be established, a presentation context for the Verification SOP Class can be negotiated, a response is received from the C-ECHO request and the association is released.

3.2.1.2 Proposed Presentation Contexts

A single presentation context will be proposed for *operator initiated communication verification* as shown in Table 2.

| Presentation Context Table | | | | | |
|----------------------------|-------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name | UID | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

Table 2: Proposed presentation context for verification

3.2.1.3 SOP Specific Conformance

Standard conformance is provided for the Verification SOP Class.

3.2.1.4 Association Termination

The association will be released upon receipt of the C-ECHO-RSP message.

3.2.2 User Initiated Image Send

3.2.2.1 Associated Real-World Activity: User Send

An operator can - via a graphical user interface - initiate sending images to a remote application entity. The associated local real-world activity is **User Send** and the remote real world activity is **Storage**. All selected images will be sent over a single association.

3.2.2.2 Proposed Presentation Contexts

Two presentation contexts will be proposed for *user initiated image send* as outlined in Table 3. Digiroebs can be configured to use Little Endian Implicit transfer syntax only.

| Presentation Context Table | | | | | | |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|--|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation | |
| Name | UID | Name | UID | | | |
| Digital X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None | |
| Digital X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | SCU | None | |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None | |

Table 3: Proposed presentation contexts for operator initiated image send

3.2.2.3 SOP Specific Conformance for DX Storage SOP Class

3.2.2.3.1 Status codes for C-STORE

The behavior when receiving C-STORE response status codes is shown in Table 4. The operator will be informed via the operator’s user interface message area.

| Status Code | Meaning | Behavior when receiving status code |
|-------------|--|--|
| | Any other status code not included in this table | An error message will be posted to the operator and an error message recorded in a log file. The images will be sent again after a configurable amount of time. All unsend images will be marked as not archived in the local database. These images will not be deleted automatically. |
| A7xx | Refused – Out of Resources | |
| A9xx | Error – Data Set does not match SOP Class | The remaining images will be sent if possible. An error message will be posted to the operator and an error message recorded in a log file. All unsend images will be marked as not archived in the local database. These images will not be deleted automatically. |
| Cxxx | Error – Cannot Understand | |
| B000 | Warning – Coercion of Data Elements | The operator will be informed after all images have been sent. The images will be marked as archived in the local database for later removal. |
| B007 | Warning – Data Set does not match SOP Class | |
| B006 | Warning – Elements Discarded | |
| 0000 | Success | |

Table 4: Behavior when receiving C-STORE response status codes (operator initiated)

Extended negotiation is not supported for the *User Send* Real-World Activity.

3.2.2.3.2 Supported IOD modules

Digiroebis uses the IOD modules shown in Table 5.

| Information Entity | Module |
|--------------------|------------------------|
| Patient | Patient |
| Study | General Study |
| Series | General Series |
| | DX Series |
| Equipment | General Equipment |
| Image | General Image |
| | Image Pixel |
| | DX Anatomy Imaged |
| | DX Image |
| | DX Detector |
| | Acquisition Context |
| | SOP Common |
| | X-Ray Acquisition Dose |
| | X-Ray Filtration |
| | DX Positioning |
| | X-Ray Generation |

Table 5: IOD modules used by Digiroebs

3.2.2.3.3 Supported DICOM Elements

3.2.2.3.3.1 Patient Module

| Attribute Name | Tag | VR | Type | Value |
|----------------------|-------------|----|------|---|
| Patient's Name | (0010,0010) | PN | 2 | Patient's full name. |
| Patient ID | (0010,0020) | LO | 2 | Primary hospital identification number or code for the patient. |
| Patient's Birth Date | (0010,0030) | DA | 2 | Birth date of the patient. |
| Patient's Sex | (0010,0040) | CS | 2 | Sex of the named patient |

3.2.2.3.3.2 General Study Module

| Attribute Name | Tag | VR | Type | Value |
|----------------------------|-------------|----|------|---|
| Study Instance UID | (0020,000D) | UI | 1 | Unique identifier for the Study. |
| Study Date | (0008,0020) | DA | 2 | Date the Study started. |
| Study Time | (0008,0030) | TM | 2 | Time the Study started. |
| Referring Physician's Name | (0008,0090) | PN | 2 | Name of the patient's referring physician. |
| Study ID | (0020,0010) | SH | 2 | User or equipment generated Study identifier. |
| Accession Number | (0008,0050) | SH | 2 | A RIS generated number that identifies the order for the Study. |
| Study Description | (0008,1030) | LO | 3 | Institution-generated description or classification of the Study (component) performed. |

3.2.2.3.3.3 General Series Module

| Attribute Name | Tag | VR | Type | Value |
|---------------------|-------------|----|------|---|
| Modality | (0008,0060) | CS | 1 | DX |
| Series Instance UID | (0020,000E) | UI | 1 | Unique identifier of the Series. |
| Series Number | (0020,0011) | IS | 2 | A number that identifies this Series. |
| Image Laterality | (0020,0062) | CS | 2C | Image Laterality R=right, L=left, U=Unpaired, B=Both |
| Series Date | (0008,0021) | DA | 3 | Date the Series started. |
| Series Time | (0008,0031) | TM | 3 | Time the Series started. |
| Series Description | (0008,103E) | LO | 3 | User provided description of the Series. |
| Body Part Examined | (0018,0015) | CS | 3 | |

3.2.2.3.3.4 DX Series Module

| Attribute Name | Tag | VR | Type | Value |
|--------------------------|-------------|----|------|------------------|
| Modality | (0008,0060) | CS | 1 | DX |
| Presentation Intent Type | (0008,0068) | CS | 1 | FOR PRESENTATION |

3.2.2.3.3.5 General Equipment Module

| Attribute Name | Tag | VR | Type | Value |
|---------------------------|-------------|----|------|--|
| Manufacturer | (0008,0070) | LO | 2 | medigration |
| Institution Name | (0008,0080) | LO | 3 | Institution where the equipment that produced the composite instances is located. Configurable. |
| Institution Address | (0008,0081) | ST | 3 | Mailing address of the institution where the equipment that produced the composite instances is located. Configurable. |
| Station Name | (0008,1010) | SH | 3 | User defined name identifying the machine that produced the composite instances. Configurable. |
| Manufacturer's Model Name | (0008,1090) | LO | 3 | Manufacturer's model name of the equipment that produced the composite instances. |
| Device Serial Number | (0018,1000) | LO | 3 | Manufacturer's serial number of the equipment that produced the composite instances. |
| Software Versions | (0018,1020) | LO | 3 | Manufacturer's designation of software version of the equipment that produced the composite instances. |
| Date of Last Calibration | (0018,1200) | DA | 3 | Date when the image acquisition device calibration was last changed. |
| Time of Last Calibration | (0018,1201) | TM | 3 | Time when the image acquisition device calibration was last changed. |

3.2.2.3.3.6 General Image Module

| Attribute Name | Tag | VR | Type | Value |
|-------------------------|-------------|----|------|---|
| Instance Number | (0020,0013) | IS | 2 | |
| Patient Orientation | (0020,0020) | CS | 2C | |
| Content Date | (0008,0023) | DA | 2C | The date the image pixel data creation started. |
| Content Time | (0008,0033) | TM | 2C | The time the image pixel data creation started. |
| Image Type | (0008,0008) | CS | 1 | Siehe Abschnitt 3.2.2.3.3.9 |
| Acquisition Date | (0008,0022) | DA | 3 | The date the acquisition of data that resulted in this image started. |
| Acquisition Time | (0008,0032) | TM | 3 | The time the acquisition of data that resulted in this image started. |
| Image Comments | (0020,4000) | LT | 3 | User-defined comments about the image. |
| Quality Control Image | (0028,0300) | CS | 3 | Indicates whether or not this image is a quality control or phantom image. YES NO |
| Burned In Annotation | (0028,0301) | CS | 1 | Siehe Abschnitt 3.2.2.3.3.9 . |
| Lossy Image Compression | (0028,2110) | CS | 1 | Siehe Abschnitt 3.2.2.3.3.9 . |

3.2.2.3.3.7 Image Pixel Module

| Attribute Name | Tag | VR | Type | Value |
|----------------------------|-------------|----|------|---------------------------------|
| Samples per Pixel | (0028,0002) | US | 1 | 1 |
| Photometric Interpretation | (0028,0004) | CS | 1 | MONOCHROME 1 |
| Rows | (0028,0010) | US | 1 | Number of rows in the image. |
| Columns | (0028,0011) | US | 1 | Number of columns in the image. |
| Bits Allocated | (0028,0100) | US | 1 | 16 |
| Bits Stored | (0028,0101) | US | 1 | 14 |
| High Bit | (0028,0102) | US | 1 | 13 |
| Pixel Representation | (0028,0103) | US | 1 | 0 |
| Pixel Data | (7FE0,0010) | OB | 1 | |
| Pixel Aspect Ratio | (0028,0034) | IS | 1C | 1\1 |

3.2.2.3.3.8 DX Anatomy Imaged Module

| Attribute Name | Tag | VR | Type | Value |
|--------------------------|-------------|----|------|---|
| Image Laterality | (0020,0062) | CS | 1 | Laterality of (possibly paired) body part examined. R = right L = left U = unpaired B = both left and right |
| Anatomic Region Sequence | (0008,2218) | SQ | 2 | |

3.2.2.3.3.9 DX Image Module

| Attribute Name | Tag | VR | Type | Value |
|-----------------------------------|-------------|----|------|-----------------------------|
| Image Type | (0008,0008) | CS | 1 | ORIGINAL\PRIMARY\ |
| Samples per Pixel | (0028,0002) | US | 1 | Siehe Abschnitt 3.2.2.3.3.6 |
| Photometric Interpretation | (0028,0004) | CS | 1 | Siehe Abschnitt 3.2.2.3.3.6 |
| Bits Allocated | (0028,0100) | US | 1 | Siehe Abschnitt 3.2.2.3.3.6 |
| Bits Stored | (0028,0101) | US | 1 | Siehe Abschnitt 3.2.2.3.3.6 |
| High Bit | (0028,0102) | US | 1 | Siehe Abschnitt 3.2.2.3.3.6 |
| Pixel Representation | (0028,0103) | US | 1 | Siehe Abschnitt 3.2.2.3.3.6 |
| Pixel Intensity Relationship | (0028,1040) | CS | 1 | LIN |
| Pixel Intensity Relationship Sign | (0028,1041) | SS | 1 | 1 |
| Rescale Intercept | (0028,1052) | DS | 1 | 0 |
| Rescale Slope | (0028,1053) | DS | 1 | 1 |
| Rescale Type | (0028,1054) | LO | 1 | US |
| Presentation LUT Shape | (2050,0020) | CS | 1 | INVERSE |
| Lossy Image Compression | (0028,2110) | CS | 1 | 00 |
| Patient Orientation | (0020,0020) | CS | 1 | |
| Calibration Image | (0050,0004) | CS | 3 | NO |
| Burned In Annotation | (0028,0301) | CS | 1 | NO |
| Window Center | (0028,1050) | DS | 1C | Window Center for display. |
| Window Width | (0028,1051) | DS | 1C | Window Width for display. |

3.2.2.3.3.10 DX Detector Module

| Attribute Name | Tag | VR | Type | Value |
|------------------------------------|-------------|----|------|--------|
| Detector Type | (0018,7004) | CS | 2 | DIRECT |
| Detector Configuration | (0018,7005) | CS | 3 | AREA |
| Detector Description | (0018,7006) | LT | 3 | |
| Detector Mode | (0018,7008) | LT | 3 | |
| Detector ID | (0018,700A) | CS | 3 | |
| Date of Last Detector Calibration | (0018,700C) | DA | 3 | |
| Time of Last Detector Calibration | (0018,700E) | TM | 3 | |
| Detector Active Time | (0018,7014) | DS | 3 | |
| Detector Conditions Nominal Flag | (0018,7000) | CS | 3 | |
| Detector Temperature | (0018,7001) | DS | 3 | |
| Imager Pixel Spacing | (0018,1164) | DS | 1 | |
| Detector Element Physical Size | (0018,7020) | DS | 3 | |
| Detector Active Shape | (0018,7024) | CS | 3 | |
| Detector Active Dimension(s) | (0018,7026) | DS | 3 | |
| Detector Manufacturer Name | (0018,702A) | CS | 3 | |
| Detector Manufacturer's Model Name | (0018,702B) | CS | 3 | |

3.2.2.3.3.11 Acquisition Context Module

| Attribute Name | Tag | VR | Type | Value |
|----------------------------|-------------|----|------|-------|
| AcquisitionContextSequence | (0040,0555) | SQ | 2 | |

3.2.2.3.3.12 SOP Common Module

| Attribute Name | Tag | VR | Type | Value |
|------------------------|-------------|----|------|---------------------------------------|
| SOP Class UID | (0008,0016) | UI | 1 | 1.2.840.10008.5.1.4.1.1.1.1 |
| SOP Instance UID | (0008,0018) | UI | 1 | Uniquely identifies the SOP Instance. |
| Specific Character Set | (0008,0005) | CS | 1C | ISO_IR 100 |

3.2.2.3.3.13 X-Ray Acquisition Dose Module

| Attribute Name | Tag | VR | Type | Value |
|-----------------------------|-------------|----|------|---|
| KVP | (0018,0060) | DS | 2 | Peak kilo voltage output of the X-Ray generator used. |
| X-Ray Tube Current | (0018,1151) | IS | 2C | X-Ray Tube Current in mA. |
| Exposure Time | (0018,1150) | IS | 2C | Duration of X-Ray exposure in msec. |
| Exposure | (0018,1152) | IS | 2C | The exposure expressed in mAs. |
| Distance Source to Detector | (0018,1110) | DS | 3 | Distance in mm from source to detector center. |
| Grid | (0018,1166) | CS | 3 | |
| Imager Pixel Spacing | (0018,1164) | DS | 1 | Siehe Abschnitt 3.2.2.3.3.10 |
| Focal Spot | (0018,1190) | DS | 3 | Nominal focal spot size in mm used to acquire this image. |
| Image Area Dose Product | (0018,115E) | DS | 3 | X-Ray dose, measured in dGy*cm ² , to which the patient was exposed for the acquisition of this image. |

3.2.2.3.3.14 X-Ray Filtration Module

| Attribute Name | Tag | VR | Type | Value |
|----------------|-------------|----|------|---|
| Filter Type | (0018,1160) | SH | 3 | STRIP WEDGE BUTTERFLY MULTIPLE NONE |

3.2.2.3.3.15 DX Positioning Module

| Attribute Name | Tag | VR | Type | Value |
|-----------------------------|-------------|----|------|--|
| View Position | (0018,5101) | CS | 3 | Radiographic view of the image relative to the imaging subject's orientation. |
| Distance Source to Detector | (0018,1110) | DS | 3 | Siehe Abschnitt 3.2.2.3.3.13 . |
| Positioner Type | (0018,1508) | CS | 2 | CARM COLUMN MAMMOGRAPHIC PANORAMIC CEPHALOSTAT RIGID NONE |
| Positioner Primary Angle | (0018,1510) | DS | 3 | Position of the X-Ray beam about the patient from the RAO to LAO direction where movement from RAO to vertical is positive, if Positioner Type (0018,1508) is CARM. |
| Detector Primary Angle | (0018,1530) | DS | 3 | Angle of the X-Ray beam in the row direction in degrees relative to the normal to the detector plane. Positive values indicate that the X-Ray beam is tilted toward higher numbered columns. Negative values indicate that the X-Ray beam is tilted toward lower numbered columns. |

Copyright © medigration GmbH. All rights reserved.
 Alle Rechte vorbehalten.

3.2.2.3.3.16 X-Ray Generation Module

| Attribute Name | Tag | VR | Type | Value |
|-----------------------------------|-------------|----|------|---------------------|
| Exposure Control Mode | (0018,7060) | CS | 3 | AUTOMATIC MANUAL |
| Exposure Control Mode Description | (0018,7062) | LT | 3 | |
| Exposure Status | (0018,7064) | CS | 3 | NORMAL ABORTED |

3.2.2.4 Association Termination

The association will be released upon receipt of the C-STORE-RSP message for the last sent image or upon receipt of refused or unknown status code.

If the peer AE aborts the association prematurely, all unsent SOP Instances are considered failed.

3.2.3 User initiated Print (Dicom Option)

3.2.3.1 Associated Real-World Activity

An operator can initiate printing images to a remote application entity. The associated local real-world activity is **User Print** and the remote real world activity is **Print**.

3.2.3.2 Proposed Presentation Contexts

A single presentation context will be proposed for *user initiated print* as shown in Table 6.

| Presentation Context Table | | | | | |
|---------------------------------------|-----------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name | UID | | |
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | SCU | None |

Table 6: Proposed presentation context for print

3.2.3.3 SOP Specific Conformance

3.2.3.3.1 Basic Grayscale Print Management

Digiroebs supports the SOP classes defined by the Basic Grayscale Print Management Meta SOP Class listed in the following table.

| SOP Class Name | SOP Class UID |
|---------------------------|------------------------|
| Basic Film Session | 1.2.840.10008.5.1.1.1 |
| Basic Film Box | 1.2.840.10008.5.1.1.2 |
| Basic Grayscale Image Box | 1.2.840.10008.5.1.1.4 |
| Printer | 1.2.840.10008.5.1.1.16 |

Table 7: SOP Classes of Basic Grayscale Print Management Meta SOP Class

The SOP specific conformance of the classes listed above is described in the subsequent sections.

3.2.3.3.2 Basic Film Session

3.2.3.3.2.1 Supported DIMSE Services

| Name | Usage | Description |
|----------|-------|---|
| N-CREATE | M/M | Creates a new film session. |
| N-SET | U/M | Updates attribute values of the film session. |
| N-DELETE | U/M | Deletes the film session including all associated film boxes, images and annotations. |
| N-ACTION | U/U | Sending of the N-ACTION will command the SCP to print all films of the film session. |

Table 8: Supported DIMSE Services for Basic Film Session

3.2.3.3.2.2 Supported SOP Class Elements

| Attribute Name | Tag | Valid Range |
|--------------------|-------------|-----------------------|
| Number of Copies | (2000,0010) | 1 - 999 |
| Print Priority | (2000,0020) | LOW MED HIGH |
| Medium Type | (2000,0030) | Up to 16 characters. |
| Film Destination | (2000,0040) | MAGAZINE PROCESSOR |
| Film Session Label | (2000,0050) | Up to 64 characters. |

Table 9: Supported SOP Class Elements for Basic Film Session

3.2.3.3.3 Basic Film Box

3.2.3.3.3.1 Supported DIMSE Services

| Name | Usage | Description |
|------|-------|-------------|
|------|-------|-------------|

| | | |
|----------|-----|--|
| N-CREATE | M/M | Creates a new film box. |
| N-SET | U/M | Updates attribute values of the film box. |
| N-DELETE | U/M | Deletes the film box including all associated images and annotations. |
| N-ACTION | U/U | Sending of the N-ACTION will command the SCP to print all films of the film box. |

Table 10: Supported DIMSE Services for Basic Film Box

3.2.3.3.2 Supported SOP Class Elements

| Attribute Name | Tag | Valid Range |
|------------------------------|-------------|---|
| Image Display Format | (2010,0010) | Details below. |
| Annotation Display Format Id | (2010,0030) | Values from configuration file. |
| Film Orientation | (2010,0040) | PORTRAIT LANDSCAPE |
| Film Size Id | (2010,0050) | Values from configuration file. Details below. |
| Magnification Type | (2010,0060) | BILINEAR ¹ |
| Smoothing Type | (2010,0080) | ignored |
| Border Density | (2010,0100) | BLACK , WHITE or an integer number (hundreds of OD). |
| Empty Image Display Density | (2010,0110) | BLACK , WHITE or an integer number (hundreds of OD). |
| Min Density | (2010,0120) | must be >= than configured min density and <= than max density and <= than configured max density |
| Max Density | (2010,0130) | must be <= than configured max density and >= than min density and >= than configured min density |
| Trim | (2010,0140) | YES, NO |
| Configuration Information | (2010,0150) | Values from configuration file ¹ . Details below. |
| Illumination | (2010,015E) | 0 – 10000 |

¹ Value can be ignored.

| | | |
|----------------------------------|-------------|---|
| Reflected Ambient Light | (2010,0160) | 0 |
| Referenced Film Session Sequence | (2010,0500) | --- |
| > Referenced SOP Class UID | (0008,1150) | 1.2.840.10008.5.1.1.1 |
| > Referenced SOP Instance UID | (0008,1155) | SOP Instance UID from creation of Basic Film Session. |
| Referenced Image Box Sequence | (2010,0510) | --- |
| > Referenced SOP Class UID | (0008,1150) | 1.2.840.10008.5.1.1.4 or 1.2.840.10008.5.1.1.4.1 |
| > Referenced SOP Instance UID | (0008,1155) | --- |

Table 11: Supported SOP Class Elements for Basic Film Box

Image Display Format

DICOM defines the Image Display Formats **STANDARD**, **ROW**, **COL**, **SLIDE**, **SUPERSLIDE** and **CUSTOM**. Digiroebs only supports the layout **STANDARD**.

3.2.3.3.4 Basic Grayscale Image Box

3.2.3.3.4.1 Supported DIMSE Services

| Name | Usage | Description |
|-------|-------|--|
| N-SET | M/M | Updates attribute values of the image box. |

Table 12: Supported DIMSE Services for Basic Grayscale Image Box

3.2.3.3.4.2 Supported SOP Class Elements

| Attribute Name | Tag | Valid Range |
|--------------------|-------------|-------------|
| Magnification Type | (2010,0060) | BILINEAR |
| Smoothing Type | (2010,0080) | ignored |

| | | |
|--------------------------------|-------------|--|
| Min Density | (2010,0120) | ignored |
| Max Density | (2010,0130) | ignored |
| Image Position | (2020,0010) | 1 - max number of images for Image Display Format |
| Polarity | (2020,0020) | NORMAL REVERSE |
| Requested Image Size | (2020,0030) | ignored |
| Basic Grayscale Image Sequence | (2020,0110) | --- |
| >Samples Per Pixel | (0028,0002) | 1 |
| >Photometric Interpretation | (0028,0004) | MONOCHROME1 MONOCHROME 1 MONOCHROME2 MONOCHROME 2 |
| >Rows | (0028,0010) | 0 < Rows < 10000 |
| >Columns | (0028,0011) | 0 < Columns < 10000 |
| >Pixel Aspect Ratio | (0028,0034) | Any pair of valid positive integers (1 to 2 ¹⁵ -1). |
| >Bits Allocated | (0028,0100) | 8 or 16 |
| >Bits Stored | (0028,0101) | 8 or 12 |
| >High Bit | (0028,0102) | BitsStored-1 |
| >Pixel Representation | (0028,0103) | 0 = unsigned |
| >Pixel Data | (7FE0,0010) | --- |

Table 13: Supported SOP Class Elements for Basic Grayscale Image Box

3.2.3.3.5 Printer

3.2.3.3.5.1 Supported DIMSE Services

| Name | Usage | Description |
|----------------|-------|--|
| N-EVENT-REPORT | M/M | The N-EVENT-REPORT is used to report the changes of the printer status in an asynchronous way. |

| | | |
|-------|-----|--|
| N-GET | U/M | The N-GET is used to retrieve information about the printer. |
|-------|-----|--|

Table 14: Supported DIMSE Services for Printer

3.2.3.3.5.2 Supported SOP Class Elements

| Attribute Name | Tag | Valid Range |
|---------------------------|-------------|------------------------------|
| Printer Status | (2110,0010) | NORMAL WARNING FAILURE |
| Printer Status Info | (2110,0020) | |
| Printer Name | (2110,0030) | |
| Manufacturer | (0008,0070) | |
| Manufacturer's Model Name | (0008,1090) | |
| Device Serial Number | (0018,1000) | |
| Software Version(s) | (0018,1020) | |
| Date of Last Calibration | (0018,1200) | |
| Time of Last Calibration | (0018,1201) | |

Table 15: Supported SOP Class Elements for Printer

3.2.4 Operator initiated query for Modality Worklist

3.2.4.1 Associated Real-World Activity

An operator can initiate sending a worklist query request to a remote application entity. The associated local real-world activity is **Query** and the remote real world activity is **Query Worklist**.

3.2.4.2 Proposed Presentation Contexts

One presentation context will be proposed for **operator initiated query for modality**

worklist as outlined in Table 16.

| Presentation Context Table | | | | | |
|--|------------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name | UID | | |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

Table 16: Proposed presentation contexts for operator initiated query for modality worklist

3.2.4.3 SOP Specific Conformance

3.2.4.3.1 Status codes for C-FIND

The behavior when receiving C-FIND response status codes is shown in Table 17.

| Status Code | Meaning | Behavior when receiving status code |
|-------------|--|--|
| | Any other status code not included in this table | The Query activity will be terminated. An error message will be recorded in a log file. |
| Cxxx | Failed operation. | |
| 0000 | Success | The Query activity will be terminated. The results will be displayed on a user interface. |
| FFxx | Pending | The Query activity will be continued. |

Table 17: Behavior when receiving C-STORE response status codes (operator initiated)

Extended negotiation is not supported for the **Query** Real-World Activity.

3.2.4.3.2 Supported DICOM Elements

The DICOM attributes sent in a C-FIND-RQ are listed in Table 18.

| Name / Module | Tag | Description |
|---------------------------------|-------------|-------------|
| SOP Common | | |
| Specific Character Set | (0008,0005) | |
| Scheduled Procedure Step | | |

| | | |
|--------------------------------------|-------------|---------------|
| Scheduled Procedure Step Sequence | (0040,0100) | |
| >Scheduled Station AE Title | (0040,0001) | |
| >Scheduled Procedure Step Start Date | (0040,0002) | Current date. |
| >Scheduled Procedure Step Start Time | (0040,0003) | |
| >Modality | (0008,0060) | DX |
| Requested Procedure | | |
| Requested Procedure Description | (0032,1060) | |
| Study Instance UID | (0020,000D) | |
| Imaging Service Request | | |
| Accession Number | (0008,0050) | |
| Referring Physician's Name | (0008,0090) | |
| Patient Identification | | |
| Patient's Name | (0010,0010) | |
| Patient ID | (0010,0020) | |
| Patient Demographic | | |
| Patient's Birth Date | (0010,0030) | |
| Patient's Sex | (0010,0040) | |

Table 18: Dicom attributes for Modality Worklist

Copyright © medigration GmbH. All rights reserved.
Alle Rechte vorbehalten.

3.2.5 Operator initiated Sending of Study Results

3.2.5.1 Associated Real-World Activity

An operator can initiate sending modality performed procedure steps to a remote application entity. The associated local real-world activity is **Send Study Result** and the remote real world activity is **Receive Study Result**.

3.2.5.2 Proposed Presentation Contexts

One presentation context will be proposed for **operator initiated sending of study results** as outlined in Table 19.

| Presentation Context Table | | | |
|----------------------------|-----------------|------|----------|
| Abstract Syntax | Transfer Syntax | Role | Extended |

| Name | UID | Name | UID | | Negotiation |
|-----------------------------------|-------------------------|---------------------------|-------------------|-----|-------------|
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

Table 19: Proposed presentation contexts for operator initiated sending of study results

3.2.5.3 SOP Specific Conformance

3.2.5.3.1 N-CREATE attributes

The DICOM attributes sent in a N-CREATE-RQ are listed in Table 20.

| Name / Module | Tag | Description |
|--|-------------|--|
| SOP Common | | |
| Specific Character Set | (0008,0005) | ISO_IR 100 |
| SOP Instance UID | (0008,0018) | |
| Performed Procedure Step Relationship | | |
| Scheduled Step Attribute Sequence | (0040,0270) | |
| >Study Instance UID | (0020,000D) | Study Instance UID of the study to be performed. |
| >Referenced Study Sequence | (0008,1110) | empty |
| >Accession Number | (0008,0050) | |
| >Requested Procedure ID | (0040,1001) | |
| >Requested Procedure Description | (0032,1060) | |
| >Scheduled Procedure Step ID | (0040,0009) | |
| >Scheduled Procedure Step Description | (0040,0007) | |
| >Scheduled Protocol Code Sequence | (0040,0008) | empty |
| Patient's Name | (0010,0010) | |
| Patient ID | (0010,0020) | |
| Patient's Birth Date | (0010,0030) | |
| Patient's Sex | (0010,0040) | |
| Referenced Patient Sequence | (0008,1120) | empty |

| Performed Procedure Step Information | | |
|---|-------------|------------------------|
| Performed Procedure Step ID | (0040,0253) | |
| Performed Station AE Title | (0040,0241) | AE Title of DigiRoebis |
| Performed Station Name | (0040,0242) | configurable |
| Performed Location | (0040,0243) | configurable |
| Performed Procedure Step Start Date | (0040,0244) | |
| Performed Procedure Step Start Time | (0040,0245) | |
| Performed Procedure Step End Date | (0040,0250) | empty |
| Performed Procedure Step End Time | (0040,0251) | empty |
| Performed Procedure Step Status | (0040,0252) | IN PROGRESS |
| Performed Procedure Step Description | (0040,0254) | |
| Performed Procedure Type Description | (0040,0255) | |
| Procedure Code Sequence | (0008,1032) | empty |
| Image Acquisition Results | | |
| Modality | (0008,0060) | DX |
| Study ID | (0020,0010) | |
| Performed Protocol Code Sequence | (0040,0260) | empty |
| Performed Series Sequence | (0040,0340) | empty |

Table 20: DICOM attributes sent in a N-CREATE-RQ

3.2.5.3.2 N-SET attributes

The DICOM attributes sent in a N-SET-RQ are listed in Table 21.

| Name / Module | Tag | Description |
|---|-------------|--------------------|
| SOP Common | | |
| SOP Instance UID | (0008,0018) | |
| Performed Procedure Step Information | | |

| | | |
|--------------------------------------|-------------|---|
| Performed Procedure Step Status | (0040,0252) | COMPLETED or DISCONTINUED |
| Performed Procedure Step Description | (0040,0254) | |
| Performed Procedure Type Description | (0040,0255) | |
| Procedure Code Sequence | (0008,1032) | empty |
| Performed Procedure Step End Date | (0040,0250) | |
| Performed Procedure Step End Time | (0040,0251) | |
| Image Acquisition Results | | |
| Performed Protocol Code Sequence | (0040,0260) | empty |
| Performed Series Sequence | (0040,0340) | Includes all the series performed. |
| >Performing Physician's Name | (0008,1050) | |
| >Protocol Name | (0018,1030) | |
| >Operator's Name | (0008,1070) | |
| >Series Instance UID | (0020,000E) | |
| >Series Description | (0008,103E) | |
| >Retrieve AE Title | (0008,0054) | |
| >Referenced Image Sequence | (0008,1140) | Includes all the images performed for this series. |
| >>Referenced SOP Class UID | (0008,1150) | |
| >>Referenced SOP Instance UID | (0008,1155) | |
| Radiation Dose | | |
| Total Time of Fluoroscopy | (0040,0300) | Total duration of X-Ray exposure during fluoroscopy in seconds (pedal time) during this Performed Procedure Step. |
| Total Number of Exposures | (0040,0301) | Total number of exposures made during this Performed Procedure Step. |
| Distance Source to Detector | (0018,1110) | Distance in mm from the source to detector center. |

| | | |
|-----------------------------|-------------|---|
| Distance Source to Entrance | (0040,0306) | Distance in mm from the source to the surface of the patient closest to the source during this Performed Procedure Step. Note: This may be an estimated value based on assumptions about the patient's body size and habitus. |
| Entrance Dose | (0040,0302) | Average entrance dose value measured in dGy at the surface of the patient during this Performed Procedure Step. Note: This may be an estimated value based on assumptions about the patient's body size and habitus. |
| Entrance Dose in mGy | (0040,8302) | Average entrance dose value measured in mGy at the surface of the patient during this Performed Procedure Step. Note: This may be an estimated value based on assumptions about the patient's body size and habitus. |
| Exposed Area | (0040,0303) | Typical dimension of the exposed area at the detector plane. If Rectangular: row dimension followed by column; if Round: diameter. Measured in mm. Note: This may be an estimated value based on assumptions about the patient's body size and habitus. |
| Image Area Dose Product | (0018,115E) | Total area-dose-product to which the patient was exposed, accumulated over the complete Performed Procedure Step and measured in dGy*cm*cm, including fluoroscopy. Notes: 1. The sum of the Image Area Dose Product of all images of a Series or a Study may not result in the actual area dose product to which the patient was exposed. 2. This may be an estimated value based on assumptions about the patient's body size and habitus. |

| | | |
|--------------------------------|-------------|--|
| Comments on Radiation Dose | (0040,0310) | User-defined comments on any special conditions related to radiation dose encountered during this Performed Procedure Step. |
| Exposure Dose Sequence | (0040,030E) | Exposure Dose Sequence will contain Total Number of Exposures (0040, 0301) items plus an item for each fluoroscopy episode not already counted as an exposure. |
| >Radiation Mode | (0018,115A) | Specifies X-Ray radiation mode. Enumerated Values: CONTINUOUS PULSED |
| >KVp | (0018,0060) | Peak kilo voltage output of the x-ray generator used. An average in the case of fluoroscopy (continuous radiation mode). |
| >X-ray Tube Current in μ A | (0018,8151) | X-ray Tube Current in μ A. An average in the case of fluoroscopy (continuous radiation mode). |
| >Exposure Time | (0018,1150) | Time of x-ray exposure or fluoroscopy in msec. |
| >Filter Type | (0018,1160) | Type of filter(s) inserted into the X-Ray beam (e.g. wedges). |
| >Filter Material | (0018,7050) | The X-Ray absorbing material used in the filter. May be multi-valued. |

Table 21: DICOM attributes sent in a N-CREATE-RQ

Extended negotiation is not supported for the **Send Study Result** Real-World Activity.

3.3 Association Acceptance Policy

Digiroebis will not accept any dicom associations.

4 Communication Profiles

4.1 Supported Communication Stacks

TCP/IP Network Communication is supported as defined in PS 3.8.

4.1.1 TCP/IP Stack

The TCP/IP stack is inherited from the underlying operating system.

4.1.1.1 Physical Network Media Support

No dependency exists on the physical network medium over which TCP/IP executes. The supported physical network media are inherited from the underlying operating system. Typical physical network media options include 10BASE-T Ethernet, 100BASE-TX Ethernet, FDDI and ATM.

5 Configuration

The following DICOM-related network parameters are configurable by the user via a graphical user interface:

- The title of the Digirolebs Application Entity.
- The AE Titles, IP Addresses and Port Numbers for all peer application entities. These parameters must be configured before associations can be initiated.
- Support by peer application entities for the Verification SOP Class. If supported, a connectivity test can be performed upon user request.
- The preferred transfer syntax for each peer application entity.

The following DICOM-related network parameters are configurable by a field service engineer for the all local activities:

- General association inactivity timeout (default 1800 seconds).
- Timeout waiting for a DIMSE request (default 1200 seconds).
- Timeout waiting for a DIMSE response (default 300 seconds).
- Maximum size of a received PDU (default 16KB).

6 Support of Extended Character Sets

The following extended character sets are supported:

ISO-IR 100 Latin Alphabet Supplementary Set No. 1 (ISO 8859-1)

Note: The DICOM default character set (ISO-IR 6) is a subset of ISO-IR 100.